

Karen M. Lyons Lab – UCLA

Adapted from Stephen G. Young Lab - UCLA

Updated on 11/22/2011

DNA Isolation from Mouse Biopsies – for genotyping

- 1) Lysis
 - a. Incubate mouse tissue in 500 μ l tail lysis buffer + 10 μ l proteinase K (10 mg/ml stock) overnight at 55°C
- 2) Removal of contaminants
 - a. 200 μ l of 8M ammonium acetate to precipitate any soluble proteins.
 - b. Incubate samples on ice for 20 minutes, invert tubes 3x, and then centrifuge at 4°C at 12kRPM for 20 minutes to pellet out the proteins.
- 3) DNA precipitation (can do this at room temperature)
 - a. Transfer 500 μ l your cold liquid lysate into 500 μ l of 100% isopropanol. Mix by repeated inversion. You should be able to see some precipitate form here. Centrifuge at 4°C for 10 minutes at high-speed. The DNA will precipitate at the bottom of the tube (you should see a small, translucent pellet).
 - i. Decant the isopropanol.
 - b. Wash the pellet in 1 ml of 70% ethanol. Centrifuge at 4°C for 10 minutes at high-speed. The DNA will precipitate at the bottom of the tube ((you should see a small, white pellet).
 - i. Decant the 70% ethanol, air-dry the pellet, and resuspend the DNA pellet in 50-100 μ l TE solution or ddH₂O.

Reagents:

- 1) Proteinase K (from Roche): stored at 10 mg/ml stocks at -20 °C
- 2) Tail lysis buffer: (10 mM Tris pH 8.0, 100 mM NaCl, 10 mM EDTA pH 8.0, 0.5% SDS)
To make 500 ml:

Stock soln	mL needed	Final conc
1 M Tris pH 8.0	5	10 mM
5 M NaCl	10	100 mM
0.5 M EDTA pH 8.0	10	10 mM
10 % SDS	25	0.5 %
dH ₂ O	450 mL	

- 3) 8M Ammonium Acetate (from Fisher Scientific #A637-500, F.W. 77.08)
To make 250 ml:
 - a) Add 154.16g of Ammonium Acetate to 150 ml of dH₂O
 - b) Stir to dissolve
 - c) Fill to 250 ml with of dH₂O
- 4) Tris-EDTA (TE) buffer (10 mM Tris pH7.5, 1 mM EDTA pH 8.0)