CARDIOVASCULAR PATHOLOGY GROSSING GUIDELINES

<u>NOTE:</u> Please review all complex or "interesting" gross specimens with Drs. Michael Fishbein, Gregory Fishbein, or Ryan Lau.

Specimen Type: PROSTHETIC CARDIAC VALVES

Procedure:

- 1. Take gross photographs and x-rays
- 2. Measure
- 3. Document presence of lesions, vegetations, calcifications
- 4. Comment on degree of pannus formation (i.e. tissue overgrowth)
- 5. For bioprosthetic valves, submit sections of the leaflets
- Log as gross only as specimens may be reviewed at a later time or kept for teaching purposes

Gross Template:

MMODAL COMMAND: PROSTHETIC HEART VALVE

It consists of a [bioprosthetic/mechanical/Have attending look at any Melody Valves or interesting heart hardware ***] valve measuring [measure in three dimensions ***] cm. [Describe type of mechanical valve, the inflow and outflow tract, fused, presence of calcifications or clot***]. The leaflets are [yellow, smooth, calcified, intact***]. Tears in the leaflets [are/are not***] identified. [Describe thickness of cusps ***] [Describe chordae tendineae, if present, noting thickening, shortening, or fusion***]. [Describe any papillary muscle that is included with the specimen***]. A gross photograph is taken. Representative sections are submitted in cassette [describe cassette submission***].

Cassette Submission: 1-2 cassettes of representative tissue (if present)

- Submit sections taken perpendicular to the surface of the valve from the point of attachment to the free margin. Submit cross sections of the chordae tendineae and the base of the papillary muscles, if present.
- Decalcify the sections if necessary.
- NOTE: Review all artificial valves with an attending pathologist.

CARDIOVASCULAR PATHOLOGY GROSSING GUIDELINES

