


## Accessing Images associated with a Case in Case Builder


This tip sheet describes how to view images that have been taken via Haiku/Canto.

When a case is built from an OpTime specimen that has a corresponding Haiku/Canto Image saved to the same encounter, the below header will appear at the top of Case Builder, which contains a hyperlink that will allow you to review the images.


1. Click the **Pathology Encounter Images** hyperlink.




2. Click the hyperlink that appears in the New Media section.



3. From here, you have two options.
  - a. You can click the chevron



- You will have to scroll within the case builder pane to see it:




or

- a. Click the hyperlink itself.



- This will open the image in its own floating window pane for review:



### Accessing Images Associated with a Case in Case Results

The images can also be found and reviewed in Case Results by navigating to the “Pathology History w/ Hx & In Process” snapshot report, which contains the same Media section that displays in Case Builder:



From here you can expand the image (click the chevron ) and scroll to review, **or** click the hyperlink text to open the image in a floating pane.