

Genitourinary Grossing Guidelines

Specimen Type: TESTIS BIOPSY

Comment: These are performed either for tumor or for analysis of spermatogenesis. They may be processed identically. It is important to treat these biopsies **GENTLY**, as testicular parenchyma is very delicate.

Gross Template:

Labeled with the patient's name (**), medical record number (**), designated "****", and received [fresh/in formalin] is a [number] [color, texture] core biopsies measuring ** cm in length x ** cm in diameter. The specimen is entirely submitted [*describe cassette submission*].

Cassette Submission: All tissue submitted