

Gynecologic Pathology Grossing Guidelines

Specimen Type: ECTOPIC/TUBAL PREGNANCY

Gross Template:

MMODAL COMMAND: "INSERT ECTOPIC TUBE"

It consist of a [*disrupted/intact****] salpingectomy measuring [***] cm in length x [***] cm in diameter. Fimbriae are [*present/absent****]. The lumen is patent and ranges from [*smallest to largest****] cm in diameter. There is a [*describe rupture- (measure length and width), note associated hemorrhage, blood clot****]. The lumen contains [*blood, embryo, chorionic villi, other- weigh and measure if appropriate****].

[*Representative sections are submitted/The specimen is entirely submitted / Document if gross photographs are taken****].

[*insert cassette summary****]

Cassette Submission: 3-4 cassettes

- Submit cross sections to demonstrate site of rupture, perpendicular sections
- Submit cross sections to demonstrate uninvolved fallopian tube
- **If no gestational sac is grossly identified, submit the entire fallopian tube sequentially**