BONE AND SOFT TISSUE PATHOLOGY GROSSING GUIDELINES

- NOTE: Document in your cassette summary which cassettes (for each specimen part) are submitted for decal
 - Ex: Representative sections are submitted (A1, A3; B2-B4; C3-C4 following decalcification).
 - Ex: Representative sections are submitted:
 - A1- bone (decal)
 - A2- skin resection margin, perpendicular
 - B2- tibia shave (decal)

Specimen Type: INTERVERTEBRAL DISC

Gross Template: Labeled with the patient's name (***), medical record number (***), designated "***", and received [fresh/in formalin] are *** of fragments of tan-white cartilaginous tissue ranging from [*** - ***cm] in maximum dimension, amounting in aggregate to *** x *** x ***cm. Representative sections are submitted in [describe cassette submission].

Cassette Submission: 1 cassette

- Submit more cassettes if lesion/mass is identified
- Decalcify if necessary