

Gynecologic Pathology Grossing Guidelines

Specimen Type: TOTAL VULVECTOMY (for tumor)

Gross Template:

MMODAL COMMAND: “INSERT VULVA”

It consists of a [*partial/total/simple/radical/subcutaneous****] vulvectomy measuring [**** x ****] cm, and is excised to a depth of [*****] cm. [*Provide orientation****]

The epidermis is remarkable for [*Describe any lesions – including size, type, borders, color, shape, distance to all margins****]. The specimen is sectioned to reveal [*describe cut surface and depth of invasion and distance of mass to deep margin****].

[*Describe remaining cut surface and presence of satellite lesion(s)****]. [*The specimen is entirely submitted/representative sections are submitted****] A gross photograph is taken.

[*insert cassette summary****]

Cassette Submission: 15-20 cassettes

- At least 3 sections of tumor, including sections showing deepest area of invasion and relationship to closest margin of resection.
- Sections of entire lateral margin of resection.
- Section of entire vaginal margin of resection.
- Sections of any skin surface with altered coloration and texture.
- Sections of uninvolved skin from right and left labia majora and minora, and clitoris.
- Sections of lymph nodes, separately labeled.
- **Take gross photograph**

