PULMONARY PATHOLOGY GROSSING GUIDELINES

Specimen Type: TRANSPLANT PNEUMONECTOMY

Gross Template:

MMODAL Command: "INSERT LUNG TRANSPLANT"

It consists of a [***] gram, [measure in three dimensions in cm***] cm [intact/previously incised***], [left/right***] pneumonectomy. The bronchus measures [***] cm in length x [***] cm in diameter.

The pleura is [smooth/glistening/roughened/remarkable for adhesions***]. There is a [minimal, moderate, extensive***] amount of anthracotic pigmentation. Sectioning reveals [describe location of emphysematous changes and fibrosis: peripheral or central; and which lobes are involved: upper, lower, etc***]. The endothelium of the vasculature is [tan-yellow and smooth or remarkable for calcifications***]. The mucosa of the bronchi is [tan-yellow, corrugated, and grossly unremarkable or remarkable fordescribe presence of lesion***].

[Describe number/size of lymph nodes if identified- indicate if hilar or intraparenchymal- and color/consistency***]. Gross photographs are taken. Representative sections are submitted.

CASSETTE SUMMARY:

- [***] bronchovascular margins, shave
- [***] upper lobe, periphery
- [***] upper lobe, central
- [***] middle lung, periphery
- [***] middle lung, central
- [***] lower lobe, periphery
- [***] lower lobe, central
- [***] representative lymph node candidates

MMODAL Command: "INSERT LUNG TRANSPLANT BILATERAL"

It consists of a [***] gram, [measure in three dimensions in cm***] cm [intact/previously incised***], [right/left***] pneumonectomy. The bronchus measures [***] cm in length x [***] cm in diameter.

The pleura is [smooth/glistening/roughened/remarkable for adhesions***]. There is a [minimal, moderate, extensive***] amount of anthracotic pigmentation. Sectioning reveals [describe location of emphysematous changes and fibrosis: peripheral or central; and which lobes are involved: upper, lower, etc***]. The endothelium of the vasculature is [tan-yellow and smooth or remarkable for calcifications***]. The mucosa of the bronchi is [tan-yellow, corrugated, and grossly unremarkable or remarkable fordescribe presence of lesion***].

Also present in the container is a [***] gram, [measure in three dimensions in***] cm [intact/previously incised***], [right/left***] pneumonectomy. The bronchus measures [***] cm in length x [***] cm in diameter.

The pleura is [smooth/glistening/roughened/remarkable for adhesions***]. There is a [minimal, moderate, extensive***] amount of anthracotic pigmentation. Sectioning reveals [describe location of emphysematous changes and fibrosis: peripheral or central; and which lobes are involved: upper, lower, etc***]. The endothelium of the vasculature is [tan-yellow and smooth or remarkable for calcifications***]. The mucosa of the bronchi is [tan-yellow, corrugated, and grossly unremarkable or remarkable fordescribe presence of lesion***].

PULMONARY PATHOLOGY GROSSING GUIDELINES

[Describe number/size of lymph nodes if identified- indicate if hilar or intraparenchymal- and color/consistency***]. Gross photographs are taken. Representative sections are submitted.

CASSETTE SUMMARY:

[***]	right lung bronchovascular margins, shave
[***]	right lung, upper lobe, periphery
[***]	right lung, upper lobe, central
[***]	right lung, middle lobe, periphery
[***]	right lung, middle lobe, central
[***]	right lung, lower lobe, periphery
[***]	right lung, lower lobe, central
[***]	right lung, representative lymph node candidates
[***]	left lung bronchovascular margins, shave
[***]	left lung, upper lobe, periphery
[***]	left lung, upper lobe, central
[***]	left lung, middle, periphery
[***]	left lung, middle, central
[***]	left lung, lower lobe, periphery
[***]	left lung, lower lobe, central
[***]	left lung, representative lymph node candidates

Cassette Submission: 8-10 cassettes per lung (16-20 total)

- Bronchovascular margins, shave
- Submit representative hilar and intraparenchymal lymph nodes
- Submit 2 sections per lobe (6 sections of each lung)
 - One section of peripheral parenchyma to include pleura
 - One section of central parenchyma
 - Submit 6 sections of left and right lungs- submit central and peripheral middle left lung parenchyma
 - Take sections to include interface between fibrosis and unremarkable lung parenchyma